

Hyperseal washers are the flexible way to feel secure

- Mould beautifully to all the body contours, whilst remaining exceptionally versatile
- Stretch, shape, mould or cut to fit around the stoma
- Suitable for all stoma types
- Treat skin with care – do not sting, even on vulnerable or sore skin
- A cost effective choice

Hyperseal hydrocolloid washers are available on prescription

Welland Hyperseal Hydrocolloid Washers				
Size	Contents	Order code	PIP code	NHS Supply Chain code
Small low profile (50mm x 2.5mm)	20	HWA 300	263-1414	GCF207
Small high profile (50mm x 4.0mm)	20	HWA 400	365-3714	GCF284
Large low profile (100mm x 2.5mm)	5	HWA 350	263-1448	GCF208



Call our free confidential Careline on **0800 036 0100** or visit www.clinimed.co.uk

CliniMed[®]
Essentials ●●●●●

Manufactured by Welland[®] a CliniMed[®] Group company. Welland products are distributed in the UK by CliniMed Ltd. Tel: 01628 850100 Fax: 01628 523579 Email: enquiries@clinimed.co.uk or visit www.clinimed.co.uk. For N. Ireland contact Fannins customer services Tel: 02890 735588 or email info@fanninni.com. CliniMed Ltd, a company registered in England number 01646927. Registered office: Cavell House, Knaves Beech Way, Loudwater, High Wycombe, Bucks HP10 9QY. Hyperseal[®] and CliniMed[®] are registered trademarks of CliniMed (Holdings) Ltd. ©2011 CliniMed Ltd. 1442/0811/1

HYPERSEAL

completely mouldable to stop even the most stubborn leaks



CliniMed[®]
Essentials ●●●●●

Hyperseal[®]
Mould it around your life

CliniMed Limited
Freepost HY241
High Wycombe
Bucks
HP10 8BR





Security that prevents leaks

Hyperseal® washers are ideal for stomas that have crevices, fistulas or scar tissue, which make the surrounding stoma site uneven. The Welland Hyperseal washer is available in low and high profile options and is easily mouldable.

Each washer is made of a consistency that allows you not only to mould, shape and stretch it to your individual needs but it also remains pliable during wear time. Even if they break, you can rejoin them as they don't lose their stickiness. Hyperseal washers are suitable for all types of stoma.

Repeated pouch changes can be harsh on skin – Hyperseal washers are gentle on vulnerable and sensitive skin, and do not sting even sore skin. By preventing pouch leaks, they can be a cost effective solution to problematic stomas.


Hyperseal washers are now available in three sizes to cater for more needs

NEW – a deeper washer with a convenient small diameter


Welland Hyperseal Hydrocolloid Washers						
Washer	Use with Convex	Small recess/scarring	Large recess	Mouldable	Stretches	Use instead of paste
Small low profile	X	✓	X	✓	✓	✓
NEW Small high profile	✓	X	✓	✓	✓	✓
Large low profile	X	X	✓	✓	✓	✓

Instructions for use


You may need more than one washer or only use part of it, depending on individual need.




1 Remove adhesive protector from washer




2 Cut or stretch washer if necessary to fit snugly around stoma




3 Fit washer over stoma



4 Mould washer to surrounding skin contours




5 Apply pouch over stoma ensuring pouch flange is adhered to washer



6 Secure pouch around stoma area



7 Make sure pouch fits securely around stoma and washer



8 The washer provides extra security against pouch leakages



Hyperseal[®]
Mould it around your life

For a free sample of:

- Hyperseal Washer Small low profile
- NEW** Hyperseal Washer Small high profile
- Hyperseal Washer Large low profile

Please complete this coupon and return it to:

CliniMed Ltd, FREEPOST HY241,
High Wycombe, Bucks. HP10 8BR
(NO STAMP REQUIRED)

You can also call our free confidential Careline on **0800 036 0100** or visit www.clinimed.co.uk

Mr/Mrs/Ms: _____ Initials: _____

Surname: _____

Address: _____

Postcode: _____

Telephone: _____

Email: _____

- I understand that this request will be handled by CliniMed Limited. I would like my details to be kept on file, so that I can be kept up to date with information about relevant new products and services.